



DiaDent Group International  
www.diadent.com | 1.877.342.3368

# Dia-ProSeal™

## Epoxy Resin-Based Root Canal Sealer




Dia-ProSeal Root Canal Sealer is a two-component epoxy resin-based system.


### ADVANTAGES

- Excellent sealing ability: low shrinkage to reduce gap between the sealer and canal wall. High adhesion to gutta percha and dentin
- Outstanding flow characteristics allow the material to reach and seal lateral canals, making it especially suitable for warm gutta percha technique
- High pH value: with the addition of calcium hydroxide, Dia-ProSeal possesses excellent antimicrobial activity
- Fast and easy preparation with no wasted material: 1:1 paste to paste manual mixing system on the mixing pad
- Low solubility: resistant to material breakdown over time with long-term dimensional stability
- Highly biocompatible: extruded material is easily resorbed to minimize periapical inflammation and sensitivity
- Excellent radiopacity and consistency
- Setting time: within 3-4 hours


### INSTRUCTIONS FOR USE




1 Press the syringe with two hands and squeeze out equal volume units of base (yellow) and catalyst (white) material in a 1:1 ratio on the mixing pad.



2 Mix well with the spatula for 10-20 seconds until they are completely mixed together and have turned into an ivory color.



3 Lightly coat the gutta percha point with the mixed material.



4 Insert the gutta percha into the canal. Take an X-ray to check the seal.

### RADIOGRAPHS



### ORDER INFORMATION

**Dia-ProSeal Regular Kit (#1003-201)**  
Contains: 1 Syringe of 16 g (6.5 g Base + 9.5 g Catalyst), 1 Mixing Pad and 1 Spatula

**Dia-ProSeal Intro Kit (#1003-202)**  
Contains: 1 Syringe of 4 g, 1 Mixing Pad and 1 Spatula

**Dia-ProSeal Spatula (#2001-6113)**  
**Dia-ProSeal Mixing Pad (#2001-6114)**