Operating Instructions  
 Gebruiksaanwijzing  
 Mode d'emploi  
 Gebrauchsanleitung  
 Brugervejledning

**DIRECTIONS FOR USE:**

GB

**APPLICATIONS:**

**1. CTT<sup>TM</sup> as a soft tissue trimmer:**

- Gingival modelling
- Freeing deep cavities
- Laceration of interradicular granulation tissue
- Papillectomy/removal of hyperplastic gingiva
- Exposure of intraosseous implant sections and retinated teeth

Any blood and pus are rinsed away, and the area is dried by a cotton sponge or by the air-syringe. The point of the CTT<sup>TM</sup> must be moved through the tissue with small "brushstroke" movements in a 30-45° degree, not forcing the point of the trimmer in the tissue or cutting too deep in only one cut.

**2. CTT<sup>TM</sup> used for opening of sulcus, e.g. for impressions:**

Any blood and pus are rinsed away, and the area is dried by a cotton sponge or by the air-syringe. The CTT<sup>TM</sup> is moved in the sulcus with an even slow sliding motion of the necessary depth. To avoid any tissue loss, it is important to preserve the epithelium of the marginal gingiva.

- CTT<sup>TM</sup> (Ceramic Tissue Trimmer) is used under rotation in the air turbine handpiece (300,000 - 500,000 RPM).
- Correctly used the kinetic energy causes a heat-coagulating effect without "burning" the tissue.
- To be used *without* any water/spray cooling.
- Using the CTT<sup>TM</sup> correctly, with all it's applications, you will find that the CTT<sup>TM</sup> is often able to replace not only electro-surgery and surgical blades, but in many cases also the use of retraction cord.

**NB:** When using CTT<sup>TM</sup> - avoid excessive contact to tooth or bone tissue.

**CLEANING:** Water and detergent.

**WARNING:** Do not use any solvent for cleaning purpose.

**STERILISATION:** Sterilise in dry heat at 180°C max. or by steam autoclave at 135°C for 30 minutes.

**MANUFACTURER:** RØNVIG Dental Mfg A/S · Gl. Vejlevej 59 · DK-8721 Daugaard · Denmark.

**CAUTION:** US Federal Law restricts this device to sale by or on the order of a dental professional.

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**HERSTELLER:** RØNVIG Dental Mfg A/S · Gl. Vejlevej 59 · DK-8721 Daugaard · Dänemark.
**PRODUCENT:** RØNVIG Dental Mfg A/S · Gl. Vejlevej 59 · DK-8721 Daugaard · Danmark.