
elite HD+



THE HIGH DEFINITION


A-Silicone for impression taking





Elite HD +, the A-Silicone which offers accurate and reliable results.

Elite HD+ is an impression material that combines the effects of **nanotechnology** with **hydrocompatibility**, thus guaranteeing high performance in a variety of clinical conditions.



Features

When hydrocompatibility becomes definition.



HYDROCOMPATIBLE

 Quick hydrocompatibility for an optimal user experience.

FLOWABILITY

▶ Nanoparticles in the fillers improve flowability in the gingival sulcus.

DETAILED

▶ The **synergistic action** between immediate hydrocompatibility and nano-technology ensures a consistent reproduction of details.

VERSATILITY

▶ A **wide range** of viscosities and setting times suitable for different impression technique.

PRACTICAL


➤ The 5:1 system for automatic mixing machines fits straight into the machine. It does not need a cartridge sleeve.

COMPLETE

A range of products that satisfies **different** clinical needs of the user.

Benefits

Hydrocompatibility and nano-technology, a synergistic action.




HIGH DEFINITION

Hydrocompatibility effect allows a quick identification of detail, therefore drastically reducing the risk of bubble formation.

FLOW PRESSURE

The **effect of the nanotecnology** also transpires when the material is put under pressure, absolutely during the impression taking Elite HD+ increases its fluidity to reach the less accessible areas.





Suggested techniques and combinations

PRODUCT	SUGGESTED TECHNIQUE	SUGGESTED COMBINATION
Elite HD+ Putty Soft Elite HD+ Maxi Putty Soft	Double step impression	Elite HD+ Regular Body Elite HD+ Light Body Elite HD+ Super Light Body
Elite HD+ Tray Material Elite HD+ Maxi Tray Material	Single step impression	Elite HD+ Regular Body Elite HD+ Light Body
Elite HD+ Monophase Elite HD+ Maxi Monophase	Monophase impression	

FOCUS



HYDROCOMPATIBILITY

The **combination of nanometric charges**, pre-treated with polar tensioactives, makes **Elite HD+** quickly hydrocompatible after the application of the fluid in the gingival sulcus.


IDEAL FLOWABILITY

The **application of nano-particles** allows **Elite HD+** to reach an ideal level of flowability within the gengival sulcus.

IMPRESSION DETAILS



Regular Body Normal Setting



Light Body Normal Setting

Codes




Elite HD+ Putty Soft - Very high viscosity A-Silicone

Code	Setting time	Packaging
C203000	Normal Set	Chan david Dank, 200 and Dane 1, 200 and Catalyint by the
C203010	Fast Set	Standard Pack: 250 ml Base + 250 ml Catalyst tubs
C203002	Normal Set	For Pools 2 of 450 or I Pools 2 of 450 or I Catalog to I had to I had
C203012	Fast Set	Eco Pack: 2 x 450 ml Base + 2 x 450 ml Catalyst tubs
C203100	Normal Set	Intro Kit: 250 ml Base + 250 ml Catalyst tubs Putty Soft Normal, 1 x 50 ml Light Normal, 1 x 50 ml Regular Normal + 12 mixing tips yellow + 6 oral tips yellow + 1 Dispenser D2



Elite HD+ Maxi Putty Soft - Very high viscosity A-Silicone






Elite HD+ Tray Material - High viscosity A-Silicone

Code	Setting time	Packaging
C202032	Fast Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 6 mixing tips green



Elite HD+ Maxi Tray Material - High viscosity A-Silicone

Code	Setting time	Packaging
C202320	Fast Set	Standard Pack: 2 x 380 ml cartridges (Base + Catalyst) + 15 dynamic-static mixing tips + 2 mixing tip lockers
C202330	Fast Set	Eco Pack: 6 x 380 ml cartridges (Base + Catalyst) + 2 mixing tip lockers



Elite HD+ Monophase - Medium viscosity A-Silicone

Code	Setting time	Packaging
C202020	Normal Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 12 mixing tips yellow

Elite HD+ Maxi Monophase - Medium viscosity A-Silicone

- 10	
5:1	

Line HD i	maxi monopii	ase medium viscosity it sincone
Code	Setting time	Packaging
C202300	Normal Set	Standard Pack: 2 x 380 ml cartridges (Base + Catalyst) + 15 dynamic-static mixing tips + 2 mixing tip lockers
C202310	Normal Set	Eco Pack: 6 x 380 ml cartridges (Base + Catalyst) + 2 mixing tip lockers

Elite HD+ Regular Body - Medium viscosity A-Silicone



Code	Setting time	Packaging
C203020	Normal Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 12 mixing tips yellow
C203025	Normal Set	Standard Pack: 2 x 90 ml tubes (Base + Catalyst)
C203100	Normal Set	Intro Kit: 250 ml Base + 250 ml Catalyst tubs Putty Soft Normal, 1 x 50 ml Light Normal, 1 x 50 ml Regular Normal + 12 mixing tips yellow + 6 intraoral tips + 1 Dispenser D2

Elite HD+ Light Body - Low viscosity A-Silicone



Code	Setting time	Packaging
C203030	Normal Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 12 mixing tips yellow
C203035	Normal Set	Standard Pack: 2 x 90 ml tubes (Base + Catalyst)
C203100	Normal Set	Intro Kit: 250 ml Base + 250 ml Catalyst tubs Putty Soft Normal, 1 x 50 ml Light Normal, 1 x 50 ml Regular Normal + 12 mixing tips yellow + 6 intraoral tips + 1 Dispenser D2

Elite HD+ Light Body - Low viscosity A-Silicone



Code	Setting time	Packaging
C203040	Fast Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 12 mixing tips yellow

Elite HD+ Super Light Body - Low viscosity A-Silicone



Code	Setting time	Packaging
C203050	Fast Set	Standard Pack: 2 x 50 ml cartridges (Base + Catalyst) + 12 mixing tips yellow





Code	Model
6000-0000	Sympress – 230 V
6000-1000	Sympress – 120 V
6000-2000	Sympress 230 V plug UK
6000-3000	Sympress 100 V

For more information about Sympress, visit our website www.zhermack.com

Accessories


























Application pictures courtesy of: Dr. M. Villaroel

Fulfilling your needs