

Cobalt Caries Detector Dye



Shelf Life:

● 2 Years From The Date Of Manufacture

✓ Components:

- 6 ml bottle * 1
- Disposable tips
- 1 ml Syringe * 1

✓ Benefits:

- Stains Tooth Hard Tissue Affected by Caries
- Remineralized Dentin Is Not Stained
- Providing Maximum Security

✓ Applications:

- Helps Minimize the Removal of Remineralized and Healthy Dentin
- Protecting the Vitality of the Pulp
- Conserving the Maximum Amount of Healthy Tissue

✓ Composition:

- PG , Acid Red



Filpex

cobalt

خمیر کلسیم هیدروکساید با یدوفرم



Shelf Life:

- 2 Years From The Date Of Manufacture

✓ Components:

- 1.1 gr syringe (1 ml) * 2
- Disposable tips

✓ Benefits:

- Anti-Bacterial Effect
- Radiopacity
- Pre-Mixed Paste
- Easy Delivery

✓ Applications:

- Apexification
- Weeping Canal
- Root Canal Filling for Deciduous Teeth
- Internal and External Root Resorption

✓ Composition:

- Calcium Hydroxide , Iodoform , Silicone Oil



Cobalt Prep



آرسی پریپ



Shelf Life:
● 2 Years From The Date Of Manufacture

✓ Components:

- 6.1 gr syringe (5 ml) * 2
- Disposable tips
- IFU

✓ Benefits:

- Removal of Inorganic Substance Through Chelation (EDTA-Ca²⁺)
- Lubricating Effect
- Easier Filling
- Produce Clean Canal
- Without Urea Peroxide
- Water Free Formulation

✓ Applications:

- Root Canal Preparation
- Presence of Calcification in the Canal
- Smear Layer Removal
- Narrow Canals
- Immediately Adhesion Is Possible

✓ Composition:

- PEG , EDTA



Aluminostat



ژل بندآورنده خون



Shelf Life:

- 2 Years From The Date Of Manufacture

✓ Components:

- 3 ml syringe * 1
- Disposable tips
- IFU

✓ Benefits:

- Designed for the Esthetic Zone
- Transparent Gel
- Leaves no Residue
- Rinses Off Easily
- Non-Drip Gel
- Does Not Interfere With Bonding

✓ Applications:

- Stops All Minor Intraoral Bleeding

✓ Composition:

- Aluminium Chloride



Cobalt Etch



ایجاد



Shelf Life:

● 2 Years From The Date Of Manufacture

✓ Components:

- 3.7gr syringe (3 ml) * 3
- Disposable tips
- IFU

✓ Benefits:

- Self-Limiting on Dentin
- With Polymeric Gelling Agent
- Rinses Completely With no Residue
- Penetrates Smallest Fissures
- Won't Run on a Vertical Surfaces
- Precise Placement

✓ Applications:

- For Etch & Rinse Technique to Etch Dentin and Enamel
- Remove the Salts Created by Etching Porcelain
- Clean the Surfaces

✓ Composition:

- 37% Phosphoric Acid



Cobalt Paste



خمیر کلیم هیدروکاید



Shelf Life:

● 2 Years From The Date Of Manufacture

✓ Components:

- 5.1 gr syringe (3 ml) * 1
- Disposable tips
- IFU

✓ Benefits:

- Water Soluble Formulation
- Rinse Off Easily
- Pre-Mixed Paste
- Optimal Radiopacity

✓ Applications:

- Apexification
- Weeping Canal
- Root Canal Filling for Deciduous Teeth
- Internal and External Root Resorption

✓ Composition:

- Calcium Hydroxide , PEG



Porcelain Etch 9%



پرسلن اچ کبالت



Porcelain Etch 6%



پرسلن اچ کبالت



Shelf Life:

● 2 Years From The Date Of Manufacture

✓ Components:

- 3.4 gr syringe (3 ml) * 1
- Disposable tips
- IFU

✓ Benefits:

- Buffered Hydrofluoric Acid Gel
- Different Concentration for Different Restorations
- Increasing the Surface Area of the Porcelain
- Microscopically Roughen the Feldspathic Porcelain or Glass Surface

✓ Applications:

- For Cementation of Ceramic Restorations
- For Intraoral Repair of Ceramics

✓ Composition:

- 5.6 % Hydrofluoric Acid (6%)
- 9.6 % Hydrofluoric Acid (9%)

